


Certificate of Analysis

General properties:

Product Name: 4 arm PEG SH
 Molecular Weight: PEG MW:750/1000/2000/3400/5000 Dalton
 Catalog Number: R-1807
 Lot Number: RJ178243

Product Structure:



Product Specifications:

Items	Specification	Results
Appearance	White/Off-white solid or viscous liquid	White/Off-white solid or viscous liquid
¹ H NMR	Identical to the structure	Identical to the structure
Substitution of SH by NMR	>90%	>90%

Date of Analysis: 2017.09.27

Recommended Retest date: 2018.09.26

Recommended Storage: -20°C

Stability: 1 year

This product is for research use only, Not for human or animal use