


Certificate of Analysis

General properties:

Product Name: pH-TP-525 荧光染料, pH:6.2
Molecular Weight: N/A
Catalog Number: R-PD-004
Lot Number: RA018164

Product Structure:



Product Specifications:

Items	Specification	Results
Appearance	Dark solid	Dark solid
¹ H NMR	Identical to the structure	Identical to the structure
Purity by HPLC	>95%	>95%
Ex/Em	407/810nm-525nm	407/810nm-525nm

Date of Analysis: 2018.08.16

Recommended Retest date: 2019.08.15

Recommended Storage: -20°C, protected from light and moisture

Stability: 1 year

This product is for research use only, Not for human or animal use