


## Certificate of Analysis

### General properties:

Product Name: PCL-Hyd-PEG-cRGD  
 Molecular Weight: PEG MW:1000/2000/3400/5000/10000/20000 Dalton  
 PCL MW:1000/2000/3000/5000/8000/10000/15000/20000/30000 Dalton  
 Catalog Number: R-PL3373  
 Lot Number: RS018201

### Product Structure:



### Product Specifications:

Items	Specification	Results
Appearance	Off-white/brown solid or viscous liquid	Off-white/brown solid or viscous liquid
<sup>1</sup> H NMR	Identical to the Structure	Identical to the Structure
PDI	<1.5	<1.5
Substitution of cRGD by NMR	>90%	>90%

Date of Analysis: 2018.09.20

Recommended Retest date: 2019.09.20

Recommended Storage: -20°C, protected from light and moisture

Stability: 1 year

**This product is for research use only, Not for human or animal use**